

Product Name: RPS8 (9W6) Rabbit Monoclonal Antibody
Catalog #: AMRe17375

Summary

Production Name	RPS8 (9W6) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% BSA.
Purification	Affinity purification

Immunogen

Gene Name	RPS8
Alternative Names	RPS8;
Gene ID	6202.0
SwissProt ID	P62241.Recombinant protein of human RPS8

Application

Dilution Ratio	WB: 1:1000
Molecular Weight	24kDa

Background


Belongs to the ribosomal protein S8e family.

Product Name: RPS8 (9W6) Rabbit Monoclonal Antibody
Catalog #: AMRe17375



Research Area

Image Data



Western blot detection of RPS8 in HeLa,CHO-K1,C6,Jurkat cell lysates using RPS8 antibody(1:1000 diluted).

Note

For research use only.