

Product Name: RUNX (7V14) Rabbit Monoclonal Antibody
Catalog #: AMRe17440

Summary

Production Name	RUNX (7V14) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% BSA.
Purification	Affinity purification

Immunogen

Gene Name	RUNX1/RUNX2/RUNX3
Alternative Names	
Gene ID	861/864/860
SwissProt ID	Q01196/Q13761/Q13950.A synthetic peptide of human RUNX1/2/3

Application


Dilution Ratio	WB: 1:2000
Molecular Weight	49kDa

Background

Product Name: RUNX (7V14) Rabbit Monoclonal Antibody
Catalog #: AMRe17440

Research Area

Image Data



Western blot detection of RUNX1/2/3 in Jurkat,Rat Brain,C6,CHO-K1 cell lysates using RUNX1/2/3 antibody(1:1000 diluted).

Note

For research use only.