

Product Name: Semaphorin 7A (9X9) Rabbit Monoclonal Antibody
Catalog #: AMRe17719

Summary

Production Name	Semaphorin 7A (9X9) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% BSA.
Purification	Affinity purification

Immunogen

Gene Name	SEMA7A
Alternative Names	
Gene ID	8482.0
SwissProt ID	O75326.Recombinant protein of human Semaphorin 7a

Application

Dilution Ratio	WB: 1:1000
Molecular Weight	75kDa


Background

Product Name: Semaphorin 7A (9X9) Rabbit Monoclonal Antibody
Catalog #: AMRe17719



Research Area

Image Data



Western blot detection of Semaphorin 7a in A549,HL-60,U251,U87-MG cell lysates using Semaphorin 7a antibody(1:1000 diluted).

Note

For research use only.