

Product Name: SKP2 (14Q16) Rabbit Monoclonal Antibody
Catalog #: AMRe17935

Summary

Production Name	SKP2 (14Q16) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
Purification	Affinity purification

Immunogen

Gene Name	SKP2
Alternative Names	FBL1; FBXL1; FLB1; p45; p45skp2; Skp2;
Gene ID	6502.0
SwissProt ID	Q13309.

Application

Dilution Ratio	WB 1:500-1:2000
Molecular Weight	48kDa

Background


Substrate recognition component of a SCF (SKP1-CUL1-F-box protein) E3 ubiquitin-protein ligase complex which mediates the ubiquitination and subsequent proteasomal degradation of target proteins involved in cell cycle progression, signal transduction and transcription. Substrate recognition component of a SCF (SKP1-CUL1-F-box protein) E3 ubiquitin-protein ligase complex which mediates the ubiquitination and subsequent proteasomal degradation of target proteins involved in cell cycle progression, signal transduction and transcription (PubMed:[11931757](http://www.uniprot.org/citations/11931757), PubMed:[12435635](http://www.uniprot.org/citations/12435635), PubMed:[12769844](http://www.uniprot.org/citations/12769844), PubMed:[12840033](http://www.uniprot.org/citations/12840033), PubMed:[15342634](http://www.uniprot.org/citations/15342634), PubMed:[15668399](http://www.uniprot.org/citations/15668399), PubMed:[15949444](http://www.uniprot.org/citations/15949444), PubMed:[16103164](http://www.uniprot.org/citations/16103164), PubMed:[16262255](http://www.uniprot.org/citations/16262255), PubMed:[16581786](http://www.uniprot.org/citations/16581786), PubMed:[16951159](http://www.uniprot.org/citations/16951159), PubMed:[17908926](http://www.uniprot.org/citations/17908926), PubMed:[17962192](http://www.uniprot.org/citations/17962192), PubMed:[22770219](http://www.uniprot.org/citations/22770219), PubMed:[32267835](http://www.uniprot.org/citations/32267835)). Specifically recognizes phosphorylated CDKN1B/p27kip and is involved in regulation of G1/S transition (By similarity). Degradation of CDKN1B/p27kip also requires CKS1. Recognizes target proteins ORC1, CDT1, RBL2, KMT2A/MLL1, CDK9, RAG2, FOXO1, UBP43, YTHDF2, and probably MYC, TOB1 and TAL1 (PubMed:[11931757](http://www.uniprot.org/citations/11931757), PubMed:[12435635](http://www.uniprot.org/citations/12435635), PubMed:[12769844](http://www.uniprot.org/citations/12769844), PubMed:[12840033](http://www.uniprot.org/citations/12840033), PubMed:[15342634](http://www.uniprot.org/citations/15342634), PubMed:[15668399](http://www.uniprot.org/citations/15668399), PubMed:[15949444](http://www.uniprot.org/citations/15949444), PubMed:[16103164](http://www.uniprot.org/citations/16103164), PubMed:[17962192](http://www.uniprot.org/citations/17962192), PubMed:[16581786](http://www.uniprot.org/citations/16581786), PubMed:[16951159](http://www.uniprot.org/citations/16951159), PubMed:[17908926](http://www.uniprot.org/citations/17908926), PubMed:[17962192](http://www.uniprot.org/citations/17962192), PubMed:[22770219](http://www.uniprot.org/citations/22770219), PubMed:[32267835](http://www.uniprot.org/citations/32267835)).

Product Name: SKP2 (14Q16) Rabbit Monoclonal Antibody
Catalog #: AMRe17935

href="http://www.uniprot.org/citations/32267835" target="_blank">32267835). Degradation of TAL1 also requires STUB1 (PubMed:17962192). Recognizes CDKN1A in association with CCNE1 or CCNE2 and CDK2 (PubMed:16262255). Promotes ubiquitination and destruction of CDH1 in a CK1-dependent manner, thereby regulating cell migration (PubMed:22770219).

Research Area

Image Data



Western blot analysis of SKP2 expression in HeLa cell lysate.

Note

For research use only.