

Product Name: SLCO5A1 (6J8) Rabbit Monoclonal Antibody
Catalog #: AMRe17964



Summary

Production Name	SLCO5A1 (6J8) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% BSA.
Purification	Affinity purification

Immunogen

Gene Name	SLCO5A1
Alternative Names	
Gene ID	81796.0
SwissProt ID	Q9H2Y9.Recombinant protein of human SLCO5A1

Application

Dilution Ratio	WB: 1:1000
Molecular Weight	92kDa


Background

Product Name: SLCO5A1 (6J8) Rabbit Monoclonal Antibody
Catalog #: AMRe17964



Research Area

Image Data



Western blot detection of SLCO5A1 in K562, HeLa cell lysates using SLCO5A1 antibody (1:1000 diluted).

Note

For research use only.