

Product Name: SOX9 (4K18) Rabbit Monoclonal Antibody
Catalog #: AMRe18143

Summary

Production Name	SOX9 (4K18) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% BSA.
Purification	Affinity purification

Immunogen

Gene Name	SOX9
Alternative Names	
Gene ID	6662.0
SwissProt ID	P48436.Recombinant protein of human SOX9

Application


Dilution Ratio	WB: 1:2000-1:10000
Molecular Weight	56kDa

Background

Product Name: SOX9 (4K18) Rabbit Monoclonal Antibody
Catalog #: AMRe18143

Research Area

Image Data



Western blot detection of SOX9 in C6,3T3,HeLa cell lysates using SOX9 antibody(1:1000 diluted).

Note

For research use only.