

Product Name: STAT5A/B (8D5) Rabbit Monoclonal Antibody
Catalog #: AMRe18364



Summary

Production Name	STAT5A/B (8D5) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% BSA.
Purification	Affinity purification

Immunogen

Gene Name	
Alternative Names	Mammary gland factor (STAT5A); Signal transducer and activator of transcription 5A; Stat 5a; STAT 5b; STAT5A; Transcription factor STAT5B;
Gene ID	
SwissProt ID	P42229/P51692 .Recombinant protein of human STAT5B

Application

Dilution Ratio	WB: 1:1000
Molecular Weight	90kDa

Product Name: STAT5A/B (8D5) Rabbit Monoclonal Antibody
Catalog #: AMRe18364




Background

Carries out a dual function: signal transduction and activation of transcription. Binds to the GAS element and activates PRL-induced transcription.

Research Area

Image Data



Western blot detection of STAT5B in C6 cell lysates using STAT5B antibody(1:1000 diluted).

Note

For research use only.