

Summary

Production Name	STC1 (4H6) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB,ELISA
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
Purification	Affinity purification

Immunogen

Gene Name	STC1
Alternative Names	Stanniocalcin; STC; STC-1; Stc1;
Gene ID	6781.0
SwissProt ID	P52823.

Application

Dilution Ratio	WB 1:500-1:2000
Molecular Weight	28kDa

Background


Product Name: STC1 (4H6) Rabbit Monoclonal Antibody
Catalog #: AMRe18377



Stimulates renal phosphate reabsorption, and could therefore prevent hypercalcemia. Stimulates renal phosphate reabsorption, and could therefore prevent hypercalcemia.

Research Area

Image Data



Western blot analysis of STC1 expression in Human Prostate lysate.

Note

For research use only.