

Product Name: TAGAP (17V10) Rabbit Monoclonal Antibody
Catalog #: AMRe18627



Summary

Production Name	TAGAP (17V10) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% BSA.
Purification	Affinity purification

Immunogen

Gene Name	TAGAP
Alternative Names	ARHGAP47; IDDM21; Tagap; TAGAP1;
Gene ID	117289.0
SwissProt ID	Q8N103.Recombinant protein of human TAGAP

Application

Dilution Ratio	WB: 1:1000
Molecular Weight	81kDa


Background

Product Name: TAGAP (17V10) Rabbit Monoclonal Antibody
Catalog #: AMRe18627

May function as a GTPase-activating protein and may play important roles during T-cell activation. May function as a GTPase-activating protein and may play important roles during T-cell activation.

Research Area

Image Data



Western blot detection of TAGAP in Rat Brain,C6,CHO-K1,HeLa cell lysates using TAGAP antibody(1:1000 diluted).

Note

For research use only.