

Product Name: TMEM43 (2M3) Rabbit Monoclonal Antibody
Catalog #: AMRe19057



Summary

Production Name	TMEM43 (2M3) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB,ELISA
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% New type preservative N and 50% glycerol. Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
Purification	Affinity purification

Immunogen

Gene Name	TMEM43
Alternative Names	ARVC5; ARVD5; EDMD7; LUMA; Tmem43; Transmembrane protein 43;
Gene ID	79188.0
SwissProt ID	Q9BTV4.

Application

Dilution Ratio	WB 1:500-1:2000
Molecular Weight	45kDa

Product Name: TMEM43 (2M3) Rabbit Monoclonal Antibody
Catalog #: AMRe19057




Background

May have an important role in maintaining nuclear envelope structure by organizing protein complexes at the inner nuclear membrane. Required for retaining emerin at the inner nuclear membrane. May have an important role in maintaining nuclear envelope structure by organizing protein complexes at the inner nuclear membrane. Required for retaining emerin at the inner nuclear membrane (By similarity).

Research Area

Image Data



Western blot analysis of TMEM43 expression in Jurkat cell lysate.

Note

For research use only.