

Product Name: YTHDC2 (5G10) Rabbit Monoclonal Antibody
Catalog #: AMRe20008



Summary

Production Name	YTHDC2 (5G10) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Supplied in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% New type preservative N and 0.05% BSA.
Purification	Affinity purification

Immunogen

Gene Name	Ythdc2
Alternative Names	
Gene ID	240255.0
SwissProt ID	B2RR83.Recombinant protein of human YTHDC2

Application

Dilution Ratio	WB: 1:1000
Molecular Weight	160kDa


Background

Product Name: YTHDC2 (5G10) Rabbit Monoclonal Antibody
Catalog #: AMRe20008



Research Area

Image Data



Western blot detection of YTHDC2 in HeLa,A549,HL-60,U251,U87-MG cell lysates using YTHDC2 antibody(1:1000 diluted).

Note

For research use only.