

## Summary

<b>Production Name</b>	GAD67 Rabbit Polyclonal Antibody
<b>Description</b>	Primary antibody
<b>Host</b>	Rabbit
<b>Application</b>	WB,IP
<b>Reactivity</b>	Human,Mouse,Rat

## Performance

<b>Conjugation</b>	Unconjugated
<b>Modification</b>	Unmodified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal Antibody
<b>Form</b>	Liquid
<b>Storage</b>	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
<b>Buffer</b>	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
<b>Purification</b>	Affinity Chromatography

## Immunogen

<b>Gene Name</b>	GAD1
<b>Alternative Names</b>	Glutamate decarboxylase 1; 67 kDa glutamic acid decarboxylase; GAD 67; GAD1; GAD; GAD67
<b>Gene ID</b>	2571
<b>SwissProt ID</b>	Q99259

## Application

<b>Dilution Ratio</b>	WB: 1/500-1/1000 IP: 1/20
<b>Molecular Weight</b>	Calculated MW: 67 kDa; Observed MW: 67 kDa

**Product Name: GAD67 Rabbit Polyclonal Antibody**  
**Catalog #: APRab00173**




## Background

Catalyzes the production of GABA.

## Research Area

Signal Transduction

## Image Data



Western blot analysis of GAD67 in (1) Jurkat lysates; (2) C6 lysates using GAD67 antibody.

## Note

For research use only.