

Product Name: STIP1 Rabbit Polyclonal Antibody
Catalog #: APRab00281



Summary

Production Name	STIP1 Rabbit Polyclonal Antibody
Description	Primary antibody
Host	Rabbit
Application	WB,IHC-P,ICC/IF,FC,IP
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification	Affinity Chromatography

Immunogen

Gene Name	STIP1
Alternative Names	Hop; P60; STI1; STI1L; STIP1; Stress induced phosphoprotein 1
Gene ID	10963
SwissProt ID	P31948

Application

Dilution Ratio	WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200 IP: 1/20 FC: 1/50-1/100
Molecular Weight	Calculated MW: 63 kDa; Observed MW: 63 kDa

Product Name: STIP1 Rabbit Polyclonal Antibody
Catalog #: APRab00281




Background

Mediates the association of the molecular chaperones HSC70 and HSP90 (HSPCA and HSPCB).

Research Area

Signal Transduction

Image Data



Western blot analysis of STIP1 in Jurkat lysates using STIP1 antibody.

Note

For research use only.