

**Product Name: PPAR gamma Rabbit Polyclonal Antibody**  
**Catalog #: APRab00391**

---

## Summary

<b>Production Name</b>	PPAR gamma Rabbit Polyclonal Antibody
<b>Description</b>	Primary antibody
<b>Host</b>	Rabbit
<b>Application</b>	WB,IHC-P
<b>Reactivity</b>	Human,Mouse,Rat

## Performance

<b>Conjugation</b>	Unconjugated
<b>Modification</b>	Unmodified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal Antibody
<b>Form</b>	Liquid
<b>Storage</b>	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
<b>Buffer</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
<b>Purification</b>	Affinity Purified

## Immunogen

<b>Gene Name</b>	PPARG CIMT1; GLM1; NR1C3; Nuclear receptor subfamily 1 group C member 3; Peroxisome proliferator activated nuclear receptor gamma variant 1; Peroxisome proliferator activated receptor gamma 1; Peroxisome Proliferator Activated Receptor gamma;
<b>Alternative Names</b>	Peroxisome proliferator-activated receptor gamma; PPAR gamma; PPAR-gamma; PPARG; PPARG_HUMAN; PPARG1; PPARG2; PPARgamma.
<b>Gene ID</b>	5468
<b>SwissProt ID</b>	P37231

## Application

<b>Dilution Ratio</b>	WB: 1/500-1/1000 IHC: 1/50-1/100
-----------------------	----------------------------------

---

**Product Name: PPAR gamma Rabbit Polyclonal Antibody**  
**Catalog #: APRab00391**

---

**Molecular Weight**

Calculated MW: 58 kDa; Observed MW: 53,57 kDa


**Background**

Receptor that binds peroxisome proliferators such as hypolipidemic drugs and fatty acids. Once activated by a ligand, the receptor binds to a promoter element in the gene for acyl-CoA oxidase and activates its transcription. It therefore controls the peroxisomal beta-oxidation pathway of fatty acids. Key regulator of adipocyte differentiation and glucose homeostasis.


**Research Area**

Epigenetics and Nuclear Signaling

**Image Data**



Immunohistochemistry analysis of paraffin-embedded mouse Heart Tissue using PPAR gamma antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Western blot analysis of PPAR gamma in Hela, 3T3, PC-12 lysates using PPAR gamma antibody

**Note**

---

**Product Name: PPAR gamma Rabbit Polyclonal  
Antibody**  
**Catalog #: APRab00391**



---

For research use only.