

Product Name: DLEU7 Rabbit Polyclonal Antibody
Catalog #: APRab00583



Summary

Production Name	DLEU7 Rabbit Polyclonal Antibody
Description	Primary antibody
Host	Rabbit
Application	IHC-P,ELISA
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Purification	Affinity Purified

Immunogen

Gene Name	DLEU7
Alternative Names	Leukemia-associated protein 7 (Deleted in lymphocytic leukemia 7)
Gene ID	220107
SwissProt ID	Q6UYE1

Application

Dilution Ratio	IHC: 1/50-1/100 ELISA: 1/10000
Molecular Weight	-

Background

-


Product Name: DLEU7 Rabbit Polyclonal Antibody
Catalog #: APRab00583




Research Area

Cell Biology

Image Data



Immunohistochemistry analysis of paraffin-embedded Human kidney using DLEU7 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunohistochemical analysis of paraffin-embedded Human tonsils using DLEU7 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

Note

For research use only.