

Product Name: SOX17 Rabbit Polyclonal Antibody
Catalog #: APRab00599



Summary

Production Name	SOX17 Rabbit Polyclonal Antibody
Description	Primary antibody
Host	Rabbit
Application	WB,ELISA
Reactivity	Human,Mouse

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Purification	Affinity Purified

Immunogen

Gene Name	SOX17
Alternative Names	SOX17; Transcription factor SOX-17
Gene ID	64321
SwissProt ID	Q9H6I2

Application

Dilution Ratio	WB: 1/500-1/1000 ELISA: 1/10000
Molecular Weight	Calculated MW: 44 kDa; Observed MW: 44 kDa

Background

Acts as transcription regulator that binds target promoter DNA and bends the DNA.


Product Name: SOX17 Rabbit Polyclonal Antibody
Catalog #: APRab00599



Research Area

Epigenetics and Nuclear Signaling

Image Data



Western blot analysis of SOX17 in Jurkat lysates using Sox17 antibody.

Note

For research use only.