

**Product Name: HspA6/A7 Rabbit Polyclonal Antibody**  
**Catalog #: APRab00608**



---

## Summary

<b>Production Name</b>	HspA6/A7 Rabbit Polyclonal Antibody
<b>Description</b>	Primary antibody
<b>Host</b>	Rabbit
<b>Application</b>	WB,ELISA
<b>Reactivity</b>	Human

## Performance

<b>Conjugation</b>	Unconjugated
<b>Modification</b>	Unmodified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal Antibody
<b>Form</b>	Liquid
<b>Storage</b>	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
<b>Buffer</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
<b>Purification</b>	Affinity Purified

## Immunogen

<b>Gene Name</b>	HSPA7
<b>Alternative Names</b>	HSPA7; HSP70B; Putative heat shock 70 kDa protein 7; Heat shock 70 kDa protein B; HSPA6; HSP70B'; Heat shock 70 kDa protein 6; Heat shock 70 kDa protein B'
<b>Gene ID</b>	3311.0
<b>SwissProt ID</b>	P48741

## Application

<b>Dilution Ratio</b>	WB: 1/500-1/1000 ELISA: 1/10000
<b>Molecular Weight</b>	Calculated MW: 40 kDa; Observed MW: 40 kDa

## Background

---

**Product Name: HspA6/A7 Rabbit Polyclonal Antibody**  
**Catalog #: APRab00608**




In cooperation with other chaperones, Hsp70s stabilize preexistent proteins against aggregation and mediate the folding of newly translated polypeptides in the cytosol as well as within organelles. These chaperones participate in all these processes through their ability to recognize nonnative conformations of other proteins.

## Research Area

Signal Transduction

## Image Data



Western blot analysis of HspA6/A7 in K562 lysates using HspA6/A7 antibody.

## Note

For research use only.