

**Product Name: Elongation Factor 1A1/2/P5 Rabbit
Polyclonal Antibody
Catalog #: APRab00625**

Summary

Production Name	Elongation Factor 1A1/2/P5 Rabbit Polyclonal Antibody
Description	Primary antibody
Host	Rabbit
Application	WB,ELISA
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Purification	Affinity Purified

Immunogen

Gene Name	EEF1A1/EEF1A2/EEF1A1P5 EEF1A1; EEF1A; EF1A; LENG7; Elongation factor 1-alpha 1; EF-1-alpha-1; Elongation factor Tu; EF-Tu; Eukaryotic elongation factor 1 A-1; eEF1A-1; Leukocyte receptor cluster member 7; EEF1A2; EEF1AL; STN; Elongation factor 1-alpha 2; EF-1-alpha-2; Eukaryotic elongation factor 1 A-2; eEF1A-2; Statin-S1; Putative elongation factor 1-alpha-like 3; EF-1-alpha-like 3; Eukaryotic elongation factor 1 A-like 3; eEF1A-like 3
Alternative Names	
Gene ID	1915/1917/158078
SwissProt ID	P68104/Q05639/Q5VTE0

Application

Dilution Ratio	WB: 1/500-1/1000 ELISA: 1/10000
-----------------------	---------------------------------

**Product Name: Elongation Factor 1A1/2/P5 Rabbit
Polyclonal Antibody**
Catalog #: APRab00625

Molecular Weight

Calculated MW: 50 kDa; Observed MW: 50 kDa


Background

This protein promotes the GTP-dependent binding of aminoacyl-tRNA to the A-site of ribosomes during protein biosynthesis.

Research Area

Epigenetics and Nuclear Signaling

Image Data



Western blot analysis of Elongation Factor 1A1/2/P5 in HepG2 lysates using Elongation Factor 1A1/2/P5 antibody.

Note

For research use only.