

Product Name: ANKRD20A1 Rabbit Polyclonal Antibody
Catalog #: APRab00637



Summary

Production Name	ANKRD20A1 Rabbit Polyclonal Antibody
Description	Primary antibody
Host	Rabbit
Application	WB,ELISA
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Purification	Affinity Purified

Immunogen

Gene Name	ANKRD20A1
Alternative Names	ANKRD20A1; ANKRD20A; Ankyrin repeat domain-containing protein 20A1
Gene ID	84210
SwissProt ID	Q5TYW2

Application

Dilution Ratio	WB: 1/500-1/1000 ELISA: 1/10000
Molecular Weight	Calculated MW: 94 kDa; Observed MW: 94 kDa

Background

Located in plasma membrane.


Product Name: ANKRD20A1 Rabbit Polyclonal Antibody
Catalog #: APRab00637



Research Area

Cell Biology

Image Data



Western blot analysis of ANKRD20A1 in 3T3 lysates using ANKRD20A1 antibody.

Note

For research use only.