

Summary

Production Name	Endoglycan Rabbit Polyclonal Antibody
Description	Primary antibody
Host	Rabbit
Application	WB,ELISA
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Purification	Affinity Purified

Immunogen

Gene Name	PODXL2
Alternative Names	PODXL2; Podocalyxin-like protein 2; Endoglycan
Gene ID	50512
SwissProt ID	Q9NZ53

Application

Dilution Ratio	WB: 1/500-1/1000 ELISA: 1/10000
Molecular Weight	Calculated MW: 65 kDa; Observed MW: 65 kDa

Background

Acts as a ligand for vascular selectins.


Product Name: Endoglycan Rabbit Polyclonal Antibody
Catalog #: APRab00640



Research Area

Immunology

Image Data



Western blot analysis of Endoglycan in HeLa lysates using Endoglycan antibody.

Note

For research use only.