

Summary

Production Name	Cripto 3 Rabbit Polyclonal Antibody
Description	Primary antibody
Host	Rabbit
Application	WB,ELISA
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Purification	Affinity Purified

Immunogen

Gene Name	TDGF1P3 TDGF1P3; CRIPTO3; TDGF2; TDGF3; Putative teratocarcinoma-derived growth factor 3;
Alternative Names	Cripto-3 growth factor; Epidermal growth factor-like cripto protein CR3; Teratocarcinoma-derived growth factor 1 pseudogene 3
Gene ID	6998.0
SwissProt ID	P51864

Application

Dilution Ratio	WB: 1/500-1/1000 ELISA: 1/10000
Molecular Weight	Calculated MW: 21 kDa; Observed MW: 21 kDa

Product Name: Cripto 3 Rabbit Polyclonal Antibody
Catalog #: APRab00641




Background

Could play a role in the determination of the epiblastic cells that subsequently give rise to the mesoderm.

Research Area

Cell Biology

Image Data



Western blot analysis of Cripto 3 in HepG2, 293T, AD293, HeLa lysates using Cripto 3 antibody.

Note

For research use only.