

Summary

Production Name	IL-22 Rabbit Polyclonal Antibody
Description	Primary antibody
Host	Rabbit
Application	WB,ELISA
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Purification	Affinity Purified

Immunogen

Gene Name	IL22
Alternative Names	IL-22; IL22; ILTIF; ZCYTO18; Interleukin-22; IL-22; Cytokine Zcyto18; IL-10-related T-cell-derived-inducible factor; IL-TIF
Gene ID	50616
SwissProt ID	Q9GZX6

Application

Dilution Ratio	WB: 1/500-1/1000 ELISA: 1/10000
Molecular Weight	Calculated MW: 20 kDa; Observed MW: 20 kDa

Background

Product Name: IL-22 Rabbit Polyclonal Antibody
Catalog #: APRab00642




Cytokine that contributes to the inflammatory response in vivo.

Research Area

Immunology

Image Data



Western blot analysis of IL-22 in mouse KIDNEY lysates using IL-22 antibody.

Note

For research use only.