

Product Name: Plasminogen Receptor Rabbit Polyclonal Antibody
Catalog #: APRab00646

Summary

Production Name	Plasminogen Receptor Rabbit Polyclonal Antibody
Description	Primary antibody
Host	Rabbit
Application	WB,ELISA
Reactivity	Human,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Purification	Affinity Purified

Immunogen

Gene Name	PLGRKT
Alternative Names	Plasminogen receptor (KT; Plg-R(KT))
Gene ID	55848
SwissProt ID	Q9HBL7

Application

Dilution Ratio	WB: 1/500-1/1000 ELISA: 1/10000
Molecular Weight	Calculated MW: 17 kDa; Observed MW: 17 kDa

Background

Product Name: Plasminogen Receptor Rabbit Polyclonal Antibody
Catalog #: APRab00646




Receptor for plasminogen.

Research Area

-

Image Data



Western blot analysis of Plasminogen Receptor in HEPG2 lysates using Plasminogen Receptor antibody.

Note

For research use only.