

Product Name: Collagen XIII alpha 1 Rabbit Polyclonal Antibody
Catalog #: APRab00808

Summary

Production Name	Collagen XIII alpha 1 Rabbit Polyclonal Antibody
Description	Primary antibody
Host	Rabbit
Application	WB,ELISA
Reactivity	Human,Mouse

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Purification	Affinity Purified

Immunogen

Gene Name	COL13A1
Alternative Names	Collagen alpha-1(XIII) chain (COLXIII A1)
Gene ID	1305
SwissProt ID	Q5TAT6

Application

Dilution Ratio	WB: 1/500-1/1000 ELISA: 1/10000
Molecular Weight	Calculated MW: 70 kDa; Observed MW: 70 kDa

Background

Product Name: Collagen XIII alpha 1 Rabbit Polyclonal Antibody
Catalog #: APRab00808




Involved in cell-matrix and cell-cell adhesion interactions that are required for normal development.

Research Area

Cell Biology

Image Data



Western blot analysis of Collagen XIII alpha 1 in HeLa, mouse kidney lysates using Collagen XIII alpha 1 antibody.

Note

For research use only.