

Product Name: Phospho-STAT3 (Ser727) Rabbit Polyclonal Antibody
Catalog #: APRab00824

Summary

Production Name	Phospho-STAT3 (Ser727) Rabbit Polyclonal Antibody
Description	Primary antibody
Host	Rabbit
Application	WB,IHC-F,IHC-P,ICC/IF,IP,ELISA
Reactivity	Human,Mouse,Rat,Monkey

Performance

Conjugation	Unconjugated
Modification	Phosphorylated
Isotype	IgG
Clonality	Polyclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Purification	Affinity Purified

Immunogen

Gene Name	STAT3
Alternative Names	STAT3; APRF; Signal transducer and activator of transcription 3; Acute-phase response factor
Gene ID	6774
SwissProt ID	P40763

Application

Dilution Ratio	WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200 IP: 1/20 ELISA: 1/10000
Molecular Weight	Calculated MW: 88 kDa; Observed MW: 85 kDa

Product Name: Phospho-STAT3 (Ser727) Rabbit Polyclonal Antibody
Catalog #: APRab00824


Background

The protein encoded by this gene is a member of the STAT protein family. In response to cytokines and growth factors, STAT family members are phosphorylated by the receptor associated kinases, and then form homo- or heterodimers that translocate to the cell nucleus where they act as transcription activators.

Research Area

Signal Transduction

Image Data



Western blot analysis of Phospho-STAT3 (Ser727) in various lysates using Phospho-STAT3 (Ser727) antibody.

Note

For research use only.