

Product Name: Annexin V Rabbit Polyclonal Antibody
Catalog #: APRab03628



Summary

Production Name	Annexin V Rabbit Polyclonal Antibody
Description	Primary antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse,Rat,Monkey

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Purification	Affinity Purified

Immunogen

Gene Name	ANXA5
Alternative Names	PP4; ANX5; ENX2; RPRGL3.
Gene ID	308
SwissProt ID	P08758

Application

Dilution Ratio	WB: 1/500-1/1000
Molecular Weight	Calculated MW: 36 kDa; Observed MW: 30 kDa

Background

This protein is an anticoagulant protein that acts as an indirect inhibitor of the thromboplastin-specific complex, which is

Product Name: Annexin V Rabbit Polyclonal Antibody
Catalog #: APRab03628




involved in the blood coagulation cascade.

Research Area

Cell Biology

Image Data



Western blot analysis of Annexin V in 3T3, HeLa, C6, COS7 and MCF-7 lysates using Annexin V antibody.

Note

For research use only.