

Product Name: GABARAP Rabbit Polyclonal Antibody
Catalog #: APRab03631



Summary

Production Name	GABARAP Rabbit Polyclonal Antibody
Description	Primary antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide, pH 7.3.
Purification	Affinity Purified

Immunogen

Gene Name	GABARAP
Alternative Names	GABARAP; GABARAPL1; GABARAPL2
Gene ID	11337
SwissProt ID	O95166

Application

Dilution Ratio	WB: 1/500-1/1000
Molecular Weight	Calculated MW: 14 kDa; Observed MW: 14,16 kDa

Background

May play a role in intracellular transport of GABA(A) receptors and its interaction with the cytoskeleton.


Product Name: GABARAP Rabbit Polyclonal Antibody
Catalog #: APRab03631




Research Area

Neuroscience

Image Data



Western blot analysis of GABARAP in A549, 3T3 and Hela lysates using GABARAP antibody.



Western blot analysis of GABARAP in rat Hippocampus, rat Cerebral Cortex, rat Brain and A549 lysates using GABARAP antibody.

Note

For research use only.