

Product Name: ABHD12B Rabbit Polyclonal Antibody
Catalog #: APRab06434



Summary

Production Name	ABHD12B Rabbit Polyclonal Antibody
Description	Rabbit Polyclonal Antibody
Host	Rabbit
Application	IHC,ELISA
Reactivity	Human,Rat,Mouse

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% New type preservative N.
Purification	Affinity purification

Immunogen

Gene Name	ABHD12B
Alternative Names	ABHD12B; C14orf29; Abhydrolase domain-containing protein 12B
Gene ID	145447.0
SwissProt ID	Q7Z5M8.The antiserum was produced against synthesized peptide derived from human ABHD12B. AA range:281-330

Application

Dilution Ratio	IHC 1:100-1:300 ELISA: 1:40000
Molecular Weight	40kD

Background


similarity:Belongs to the serine esterase family.,similarity:Belongs to the serine esterase family.,

Product Name: ABHD12B Rabbit Polyclonal Antibody
Catalog #: APRab06434



Research Area

Image Data



Immunohistochemistry analysis of paraffin-embedded human colon carcinoma, using ABHD12B Antibody. The picture on the right is blocked with the synthesized peptide.

Note

For research use only.