

Product Name: ABHD14A Rabbit Polyclonal Antibody
Catalog #: APRab06435



Summary

Production Name	ABHD14A Rabbit Polyclonal Antibody
Description	Rabbit Polyclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% New type preservative N.
Purification	Affinity purification

Immunogen

Gene Name	ABHD14A
Alternative Names	ABHD14A; Abhydrolase domain-containing protein 14A
Gene ID	25864.0
SwissProt ID	Q9BUJ0.The antiserum was produced against synthesized peptide derived from human ABHD14A. AA range:201-250

Application

Dilution Ratio	WB 1:500-1:2000. ELISA: 1:40000.
Molecular Weight	30kD

Background

caution:It is uncertain whether Met-1 or Met-30 is the initiator.,function:Possible role in granule neuron


Product Name: ABHD14A Rabbit Polyclonal Antibody
Catalog #: APRab06435



development.,similarity:Belongs to the AB hydrolase superfamily. ABHD14 family.,caution:It is uncertain whether Met-1 or Met-30 is the initiator.,function:Possible role in granule neuron development.,similarity:Belongs to the AB hydrolase superfamily. ABHD14 family.,

Research Area

Image Data



Western blot analysis of lysates from RAW264.7 cells, using ABHD14A Antibody. The lane on the right is blocked with the synthesized peptide.

Note

For research use only.