

Product Name: ABHDD Rabbit Polyclonal Antibody
Catalog #: APRab06446



Summary

Production Name	ABHDD Rabbit Polyclonal Antibody
Description	Rabbit Polyclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, and 0.02% New type preservative N.
Purification	Affinity purification

Immunogen

Gene Name	ABHD13 C13orf6
Alternative Names	
Gene ID	84945.0
SwissProt ID	Q7L211.Synthesized peptide derived from human protein . at AA range: 120-200

Application

Dilution Ratio	WB 1:500-2000 ELISA 1:5000-20000
Molecular Weight	37kD

Background


sequence caution:Translated as Leu.,similarity:Belongs to the serine esterase family.,sequence caution:Translated as Leu.,similarity:Belongs to the serine esterase family.,

Product Name: ABHDD Rabbit Polyclonal Antibody
Catalog #: APRab06446



Research Area

Image Data



Western blot analysis of lysates from SH-SY5Y cells, primary antibody was diluted at 1:1000, 4°over night

Note

For research use only.