

Product Name: ALDH3B1 Rabbit Polyclonal Antibody
Catalog #: APRab06765



Summary

Production Name	ALDH3B1 Rabbit Polyclonal Antibody
Description	Rabbit Polyclonal Antibody
Host	Rabbit
Application	IHC,WB,
Reactivity	Human,Rat,Mouse

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% New type preservative N.
Purification	Affinity purification

Immunogen

Gene Name	ALDH3B1
Alternative Names	ALDH3B1; ALDH7; Aldehyde dehydrogenase family 3 member B1; Aldehyde dehydrogenase 7
Gene ID	221.0
SwissProt ID	P43353.The antiserum was produced against synthesized peptide derived from human ALDH3B1. AA range:51-100

Application

Dilution Ratio	WB 1:500 - 1:2000. IHC 1:100 - 1:300. ELISA: 1:40000..
Molecular Weight	52kD

Product Name: ALDH3B1 Rabbit Polyclonal Antibody
Catalog #: APRab06765




Background

This gene encodes a member of the aldehyde dehydrogenase protein family. Aldehyde dehydrogenases are a family of isozymes that may play a major role in the detoxification of aldehydes generated by alcohol metabolism and lipid peroxidation. The encoded protein is able to oxidize long-chain fatty aldehydes in vitro, and may play a role in protection from oxidative stress. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Feb 2014], catalytic activity: An aldehyde + NAD(P)(+) + H₂O = an acid + NAD(P)H., pathway: Alcohol metabolism; ethanol degradation; acetate from ethanol: step 2/2., similarity: Belongs to the aldehyde dehydrogenase family., tissue specificity: Highest expression in kidney and lung.,


Research Area

Glycolysis / Gluconeogenesis; Histidine metabolism; Tyrosine metabolism; Phenylalanine metabolism; Metabolism of xenobiotics by cytochrome P450; Drug metabolism;

Image Data




Immunohistochemistry analysis of paraffin-embedded human breast carcinoma tissue, using ALDH3B1 Antibody. The picture on the right is blocked with the synthesized peptide.




Western blot analysis of lysates from HUVEC cells, using ALDH3B1 Antibody. The lane on the right is blocked with the synthesized peptide.

Product Name: ALDH3B1 Rabbit Polyclonal Antibody
Catalog #: APRab06765



Western blot analysis of the lysates from RAW264.7 cells using ALDH3B1 antibody.



Western Blot analysis of various cells using ALDH3B1 Polyclonal Antibody diluted at 1: 500

Note

For research use only.