

Product Name: AVP Receptor V3 Rabbit Polyclonal Antibody
Catalog #: APRab07380

Summary

Production Name	AVP Receptor V3 Rabbit Polyclonal Antibody
Description	Rabbit Polyclonal Antibody
Host	Rabbit
Application	WB,ELISA
Reactivity	Human,Rat,Mouse

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% New type preservative N.
Purification	Affinity purification

Immunogen

Gene Name	AVPR1B
Alternative Names	AVPR1B; AVPR3; VPR3; Vasopressin V1b receptor; V1bR; AVPR V1b; AVPR V3; Antidiuretic hormone receptor 1b; Vasopressin V3 receptor
Gene ID	553.0
SwissProt ID	P47901.The antiserum was produced against synthesized peptide derived from human AVPR1B. AA range:275-324

Application

Dilution Ratio	WB 1:500 - 1:2000. ELISA: 1:20000. Not yet tested in other applications.
Molecular Weight	47kD

Product Name: AVP Receptor V3 Rabbit Polyclonal Antibody
Catalog #: APRab07380




Background

The protein encoded by this gene acts as receptor for arginine vasopressin. This receptor belongs to the subfamily of G-protein coupled receptors which includes AVPR1A, V2R and OXT receptors. Its activity is mediated by G proteins which stimulate a phosphatidylinositol-calcium second messenger system. The receptor is primarily located in the anterior pituitary, where it stimulates ACTH release. It is expressed at high levels in ACTH-secreting pituitary adenomas as well as in bronchial carcinoids responsible for the ectopic ACTH syndrome. A spliced antisense transcript of this gene has been reported but its function is not known. [provided by RefSeq, Jul 2008],function:Receptor for arginine vasopressin. The activity of this receptor is mediated by G proteins which activate a phosphatidyl-inositol-calcium second messenger system.,similarity:Belongs to the G-protein coupled receptor 1 family.,


Research Area

Calcium;Neuroactive ligand-receptor interaction;Vascular smooth muscle contraction;

Image Data




Western blot analysis of lysates from COLO and HepG2 cells, using AVPR1B Antibody. The lane on the right is blocked with the synthesized peptide.



Western blot analysis of the lysates from HUVEC cells using AVPR1B antibody.

Product Name: AVP Receptor V3 Rabbit Polyclonal Antibody
Catalog #: APRab07380



Western Blot analysis of various cells using AVP Receptor V3 Polyclonal Antibody

Note

For research use only.