

Product Name: C1QL4 Rabbit Polyclonal Antibody
Catalog #: APRab07728



Summary

Production Name	C1QL4 Rabbit Polyclonal Antibody
Description	Rabbit Polyclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, and 0.02% New type preservative N.
Purification	Affinity purification

Immunogen

Gene Name	C1QL4
Alternative Names	
Gene ID	338761.0
SwissProt ID	Q86Z23.Synthesized peptide derived from part region of human protein AA range: 110-150

Application

Dilution Ratio	WB 1:500-2000 ELISA 1:5000-20000
Molecular Weight	26kD

Background

similarity:Contains 1 C1q domain.,similarity:Contains 1 collagen-like domain.,similarity:Contains 1 C1q


Product Name: C1QL4 Rabbit Polyclonal Antibody
Catalog #: APRab07728



domain.,similarity:Contains 1 collagen-like domain.,

Research Area

Image Data



Western blot analysis of lysates from PC12 cells, primary antibody was diluted at 1:1000, 4° over night

Note

For research use only.