

Product Name: CY561 Rabbit Polyclonal Antibody
Catalog #: APRab09569



Summary

Production Name	CY561 Rabbit Polyclonal Antibody
Description	Rabbit Polyclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, and 0.02% New type preservative N.
Purification	Affinity purification

Immunogen

Gene Name	CYB561
Alternative Names	
Gene ID	1534.0
SwissProt ID	P49447.Synthesized peptide derived from human protein . at AA range: 30-110

Application

Dilution Ratio	WB 1:500-2000 ELISA 1:5000-20000
Molecular Weight	27kD

Background

cofactor: Binds 2 heme groups non-covalently.,function: Secretory vesicle-specific electron transport protein.,similarity: Contains 1 cytochrome b561 domain.,cofactor: Binds 2 heme groups non-covalently.,function: Secretory


Product Name: CY561 Rabbit Polyclonal Antibody
Catalog #: APRab09569



vesicle-specific electron transport protein.,similarity:Contains 1 cytochrome b561 domain.,

Research Area

Image Data



Western blot analysis of lysates from U2OS cells, primary antibody was diluted at 1:1000, 4° over night

Note

For research use only.