
Summary

Production Name	DACA-1 Rabbit Polyclonal Antibody
Description	Rabbit Polyclonal Antibody
Host	Rabbit
Application	WB,ELISA
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% New type preservative N.
Purification	Affinity purification

Immunogen

Gene Name	YTHDF1
Alternative Names	YTHDF1; C20orf21; YTH domain family protein 1; Dermatomyositis associated with cancer putative autoantigen 1; DACA-1
Gene ID	54915.0
SwissProt ID	Q9BYJ9.The antiserum was produced against synthesized peptide derived from human YTHDF1. AA range:1-50

Application

Dilution Ratio	WB 1:500 - 1:2000. ELISA: 1:5000
Molecular Weight	61kD

Product Name: DACA-1 Rabbit Polyclonal Antibody
Catalog #: APRab09776




Background


similarity:Contains 1 YTH domain.,similarity:Contains 1 YTH domain.,

Research Area

Image Data



Western blot analysis of lysates from 293 cells, using YTHDF1 Antibody. The lane on the right is blocked with the synthesized peptide.



Western Blot analysis of various cells using DACA-1 Polyclonal Antibody

Note

For research use only.