

Summary

Production Name	eIF2B δ Rabbit Polyclonal Antibody
Description	Rabbit Polyclonal Antibody
Host	Rabbit
Application	WB,ELISA
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% New type preservative N.
Purification	Affinity purification

Immunogen

Gene Name	EIF2B4
Alternative Names	EIF2B4; EIF2BD; Translation initiation factor eIF-2B subunit delta; eIF-2B GDP-GTP exchange factor subunit delta
Gene ID	8890.0
SwissProt ID	Q9UI10.The antiserum was produced against synthesized peptide derived from human EIF2B4. AA range:226-275

Application


Dilution Ratio	WB 1:500 - 1:2000. ELISA: 1:5000.
Molecular Weight	57kD

Background

Eukaryotic initiation factor 2B (EIF2B), which is necessary for protein synthesis, is a GTP exchange factor composed of five different subunits. The protein encoded by this gene is the fourth, or delta, subunit. Defects in this gene are a cause of leukoencephalopathy with vanishing white matter (VWM) and ovarioleukodystrophy. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2008],disease:Defects in EIF2B4 are a cause of leukodystrophy with vanishing white matter (VWM) [MIM:603896]. VWM is a leukodystrophy that occurs mainly in children. Neurological signs include progressive cerebellar ataxia, spasticity, inconstant optic atrophy and relatively preserved mental abilities. The disease is chronic-progressive with, in most individuals, additional episodes of rapid deterioration following febrile infections or minor head trauma. While childhood onset is the most common form of the disorder, some severe forms are apparent at birth. A severe, early-onset form seen among the Cree and Chippewyan populations of Quebec and Manitoba is called Cree leukoencephalopathy. Milder forms may not become evident until adolescence or adulthood. Some females with milder forms of the disease who survive to adolescence exhibit ovarian dysfunction. This variant of the disorder is called ovarioleukodystrophy.,function:Catalyzes the exchange of eukaryotic initiation factor 2-bound GDP for GTP.,similarity:Belongs to the EIF-2B alpha/beta/delta subunits family.,subunit:Complex of five different subunits; alpha, beta, gamma, delta and epsilon.,

Research Area

Image Data



Western blot analysis of EIF2B4 Antibody. The lane on the right is blocked with the EIF2B4 peptide.

Product Name: eIF2B δ Rabbit Polyclonal Antibody
Catalog #: APRab10367



Western blot analysis of the lysates from HeLa cells using EIF2B4 antibody.

Note

For research use only.