

Product Name: EPG5 Rabbit Polyclonal Antibody
Catalog #: APRab10509



Summary

Production Name	EPG5 Rabbit Polyclonal Antibody
Description	Rabbit Polyclonal Antibody
Host	Rabbit
Application	IHC
Reactivity	Human,Rat,Mouse

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% New type preservative N.
Purification	Affinity purification

Immunogen

Gene Name	EPG5 KIAA1632
Alternative Names	Ectopic P granules protein 5 homolog
Gene ID	57724.0
SwissProt ID	Q9HCE0.Synthetic peptide from human protein at AA range: 141-190

Application

Dilution Ratio	IHC 1:200-500
Molecular Weight	

Background

This gene encodes a large coiled coil domain-containing protein that functions in autophagy during starvation conditions.


Product Name: EPG5 Rabbit Polyclonal Antibody
Catalog #: APRab10509




Mutations in this gene cause Vici syndrome. [provided by RefSeq, Aug 2015],sequence caution:Translated as Trp.,similarity:Belongs to the UPF0493 family.,

Research Area

Image Data



Immunohistochemical analysis of paraffin-embedded human-liver, antibody was diluted at 1:200



Immunohistochemical analysis of paraffin-embedded human-liver, antibody was diluted at 1:200

Note

For research use only.