

Product Name: FAM3D Rabbit Polyclonal Antibody
Catalog #: APRab10815



Summary

Production Name	FAM3D Rabbit Polyclonal Antibody
Description	Rabbit Polyclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human,Rat,Mouse

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% New type preservative N.
Purification	Affinity purification

Immunogen

Gene Name	FAM3D
Alternative Names	Protein FAM3D
Gene ID	131177.0
SwissProt ID	Q96BQ1.The antiserum was produced against synthesized peptide derived from the Internal region of human FAM3D. AA range:121-170

Application

Dilution Ratio	WB 1:500-2000, ELISA 1:10000-20000
Molecular Weight	22kD

Background

similarity:Belongs to the FAM3 family.,tissue specificity:Abundantly expressed in placenta and weakly expressed in small


Product Name: FAM3D Rabbit Polyclonal Antibody
Catalog #: APRab10815



intestine.,similarity:Belongs to the FAM3 family.,tissue specificity:Abundantly expressed in placenta and weakly expressed in small intestine.,

Research Area

Image Data



Western Blot analysis of HeLa cells using FAM3D Polyclonal Antibody diluted at 1:500. Secondary antibody was diluted at 1:20000

Note

For research use only.