

Product Name: IFT1B Rabbit Polyclonal Antibody
Catalog #: APRab12410



Summary

Production Name	IFT1B Rabbit Polyclonal Antibody
Description	Rabbit Polyclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human,Rat,Mouse

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, and 0.02% New type preservative N.
Purification	Affinity purification

Immunogen

Gene Name	IFIT1B IFIT1L
Alternative Names	
Gene ID	439996.0
SwissProt ID	Q5T764.Synthesized peptide derived from human protein . at AA range: 10-90

Application

Dilution Ratio	WB 1:500-2000 ELISA 1:5000-20000
Molecular Weight	52kD

Background


similarity:Belongs to the IFIT family.,similarity:Contains 8 TPR repeats.,similarity:Belongs to the IFIT family.,similarity:Contains 8 TPR repeats.,

Product Name: IFT1B Rabbit Polyclonal Antibody
Catalog #: APRab12410



Research Area

Image Data



Western Blot analysis of Hela lysis, using primary antibody at 1:1000 dilution. Secondary antibody was diluted at 1:10000

Note

For research use only.