

Summary

Production Name	M6A Rabbit Polyclonal Antibody
Description	Rabbit Polyclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% New type preservative N.
Purification	Affinity purification

Immunogen

Gene Name	GPM6A
Alternative Names	GPM6A; M6A; Neuronal membrane glycoprotein M6-a; M6a
Gene ID	2823.0
SwissProt ID	P51674.Synthesized peptide derived from M6A . at AA range: 140-220

Application

Dilution Ratio	WB 1:500-1:2000. ELISA: 1:5000.
Molecular Weight	31kD

Background


similarity:Belongs to the myelin proteolipid protein family.,similarity:Belongs to the myelin proteolipid protein family.,

Product Name: M6A Rabbit Polyclonal Antibody
Catalog #: APRab13538



Research Area

Image Data



Western Blot analysis of extracts from rat stomach, using M6A Polyclonal Antibody.. Secondary antibody was diluted at 1:20000

Note

For research use only.