

Product Name: MPP9 Rabbit Polyclonal Antibody
Catalog #: APRab14064



Summary

Production Name	MPP9 Rabbit Polyclonal Antibody
Description	Rabbit Polyclonal Antibody
Host	Rabbit
Application	WB,
Reactivity	Human,Rat,Mouse

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% New type preservative N.
Purification	Affinity purification

Immunogen

Gene Name	MPHOSPH9
Alternative Names	MPHOSPH9; MPP9; M-phase phosphoprotein 9
Gene ID	10198.0
SwissProt ID	Q99550.The antiserum was produced against synthesized peptide derived from human MPHOSPH9. AA range:539-588

Application

Dilution Ratio	WB 1:500-2000
Molecular Weight	78kD

Background


Product Name: MPP9 Rabbit Polyclonal Antibody
Catalog #: APRab14064




PTM:Phosphorylated in M (mitotic) phase.,PTM:Phosphorylated in M (mitotic) phase.,

Research Area


Image Data



Western blot analysis of MPHOSPH9 Antibody. The lane on the right is blocked with the MPHOSPH9 peptide.



Western blot analysis of the lysates from COLO205 cells using MPHOSPH9 antibody.



Western Blot analysis of various cells using MPP9 Polyclonal Antibody

Product Name: MPP9 Rabbit Polyclonal Antibody
Catalog #: APRab14064



Note

For research use only.