

Product Name: NLE1 Rabbit Polyclonal Antibody
Catalog #: APRab14740



Summary

Production Name	NLE1 Rabbit Polyclonal Antibody
Description	Rabbit Polyclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human,Monkey

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% New type preservative N.
Purification	Affinity purification

Immunogen

Gene Name	NLE1
Alternative Names	NLE1; HUSSY-07; Notchless protein homolog 1
Gene ID	54475.0
SwissProt ID	Q9NVX2.The antiserum was produced against synthesized peptide derived from human NLE1. AA range:301-350

Application

Dilution Ratio	WB 1:500-1:2000. ELISA: 1:40000.
Molecular Weight	48kD

Background


similarity:Contains 8 WD repeats.,similarity:Contains 8 WD repeats.,

Product Name: NLE1 Rabbit Polyclonal Antibody
Catalog #: APRab14740



Research Area

Image Data



Western blot analysis of lysates from COS cells, using NLE1 Antibody. The lane on the right is blocked with the synthesized peptide.

Note

For research use only.