
Summary

Production Name	OB-cadherin Rabbit Polyclonal Antibody
Description	Rabbit Polyclonal Antibody
Host	Rabbit
Application	WB,ELISA
Reactivity	Human,Mouse

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% New type preservative N.
Purification	Affinity purification

Immunogen

Gene Name	CDH11
Alternative Names	CDH11; Cadherin-11; OSF-4; Osteoblast cadherin; OB-cadherin
Gene ID	1009.0
SwissProt ID	P55287.The antiserum was produced against synthesized peptide derived from human CDH11. AA range:421-470

Application

Dilution Ratio	WB 1:500 - 1:2000. ELISA: 1:20000. Not yet tested in other applications.
Molecular Weight	85kD

Background

Product Name: OB-cadherin Rabbit Polyclonal Antibody
Catalog #: APRab15087




This gene encodes a type II classical cadherin from the cadherin superfamily, integral membrane proteins that mediate calcium-dependent cell-cell adhesion. Mature cadherin proteins are composed of a large N-terminal extracellular domain, a single membrane-spanning domain, and a small, highly conserved C-terminal cytoplasmic domain. Type II (atypical) cadherins are defined based on their lack of a HAV cell adhesion recognition sequence specific to type I cadherins. Expression of this particular cadherin in osteoblastic cell lines, and its upregulation during differentiation, suggests a specific function in bone development and maintenance. [provided by RefSeq, Jul 2008],function: Cadherins are calcium dependent cell adhesion proteins.,function: Cadherins are calcium dependent cell adhesion proteins. They preferentially interact with themselves in a homophilic manner in connecting cells; cadherins may thus contribute to the sorting of heterogeneous cell types.,similarity: Contains 5 cadherin domains.,tissue specificity: Expressed mainly in brain but also found in other tissues. Expressed in neuroblasts.,


Research Area

Adherens_Junction

Image Data




Western blot analysis of lysates from Jurkat cells, using CDH11 Antibody. The lane on the right is blocked with the synthesized peptide.



Western blot analysis of the lysates from HeLa cells using Cytochrome P450 1A1/2 antibody.

Product Name: OB-cadherin Rabbit Polyclonal Antibody
Catalog #: APRab15087



Western Blot analysis of Jurkat cells using OB-cadherin Polyclonal Antibody

Note

For research use only.