

Product Name: Olfactory receptor 4C3 Rabbit Polyclonal Antibody
Catalog #: APRab15221



Summary

Production Name	Olfactory receptor 4C3 Rabbit Polyclonal Antibody
Description	Rabbit Polyclonal Antibody
Host	Rabbit
Application	IF, WB, ELISA
Reactivity	Human, Rat, Mouse

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% New type preservative N.
Purification	Affinity purification

Immunogen

Gene Name	OR4C3
Alternative Names	OR4C3; Olfactory receptor 4C3; Olfactory receptor OR11-98
Gene ID	256144.0
SwissProt ID	Q8NH37. The antiserum was produced against synthesized peptide derived from human OR4C3. AA range: 196-245

Application

Dilution Ratio	WB 1:500 - 1:2000. IF 1:200 - 1:1000. ELISA: 1:20000. Not yet tested in other applications.
Molecular Weight	36kD

Product Name: Olfactory receptor 4C3 Rabbit Polyclonal Antibody
Catalog #: APRab15221




Background

olfactory receptor family 4 subfamily C member 3(OR4C3) Homo sapiens Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms. [provided by RefSeq, Jul 2008],function:Odorant receptor ,similarity:Belongs to the G-protein coupled receptor 1 family,.


Research Area

Olfactory transduction;

Image Data




Immunofluorescence analysis of A549 cells, using OR4C3 Antibody. The picture on the right is blocked with the synthesized



Western blot analysis of lysates from K562 cells, using OR4C3 Antibody. The lane on the right is blocked with the synthesized

**Product Name: Olfactory receptor 4C3 Rabbit Polyclonal
Antibody**
Catalog #: APRab15221



Western blot analysis of the lysates from COLO205 cells using OR4C3 antibody.

Note

For research use only.