

Product Name: Olfactory receptor 8H3 Rabbit Polyclonal Antibody
Catalog #: APRab15329

Summary

| | |
|------------------------|---|
| Production Name | Olfactory receptor 8H3 Rabbit Polyclonal Antibody |
| Description | Rabbit Polyclonal Antibody |
| Host | Rabbit |
| Application | IF,ELISA |
| Reactivity | Human,Rat,Mouse |

Performance

| | |
|---------------------|--|
| Conjugation | Unconjugated |
| Modification | Unmodified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Form | Liquid |
| Storage | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles. |
| Buffer | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% New type preservative N. |
| Purification | Affinity purification |

Immunogen

| | |
|--------------------------|--|
| Gene Name | OR8H3 |
| Alternative Names | OR8H3; Olfactory receptor 8H3; Olfactory receptor OR11-172 |
| Gene ID | 390152.0 |
| SwissProt ID | Q8N146.The antiserum was produced against synthesized peptide derived from human OR8H3. AA range:262-311 |

Application

| | |
|-------------------------|----------------------------------|
| Dilution Ratio | IF 1:200-1:1000. ELISA: 1:10000. |
| Molecular Weight | |

Background


Product Name: Olfactory receptor 8H3 Rabbit Polyclonal Antibody
Catalog #: APRab15329

olfactory receptor family 8 subfamily H member 3(OR8H3) Homo sapiens Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms. [provided by RefSeq, Jul 2008],function:Odorant receptor ,similarity:Belongs to the G-protein coupled receptor 1 family.,

Research Area

Olfactory transduction;

Image Data



Immunofluorescence analysis of HUVEC cells, using OR8H3 Antibody. The picture on the right is blocked with the synthesized peptide.

Note

For research use only.