

Product Name: OSR2 Rabbit Polyclonal Antibody
Catalog #: APRab15522



Summary

Production Name	OSR2 Rabbit Polyclonal Antibody
Description	Rabbit Polyclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% New type preservative N.
Purification	Affinity purification

Immunogen

Gene Name	OSR2
Alternative Names	OSR2; Protein odd-skipped-related 2
Gene ID	116039.0
SwissProt ID	Q8N2R0.Synthesized peptide derived from OSR2 . at AA range: 60-140

Application

Dilution Ratio	WB 1:500-1:2000. ELISA: 1:40000.
Molecular Weight	35kD

Background

odd-skipped related transcription factor 2(OSR2) Homo sapiens OSR2 is a mammalian homolog of the Drosophila odd-skipped family of transcription factors (Lan et al., 2004 [PubMed 15175245]).[supplied by OMIM, Mar


Product Name: OSR2 Rabbit Polyclonal Antibody
Catalog #: APRab15522



2008],similarity:Belongs to the Odd C2H2-type zinc-finger protein family.,similarity:Contains 5 C2H2-type zinc fingers.,

Research Area

Image Data



Western Blot analysis of JK cells using OSR2 Polyclonal Antibody cells nucleus extracted by Minute TM Cytoplasmic and Nuclear Fractionation kit (SC-003, Invent biotech, MN, USA) .

Note

For research use only.