

**Product Name: PI 3-Kinase p85 $\beta$  Rabbit Polyclonal Antibody**  
**Catalog #: APRab16105**



---

## Summary

<b>Production Name</b>	PI 3-Kinase p85 $\beta$ Rabbit Polyclonal Antibody
<b>Description</b>	Rabbit Polyclonal Antibody
<b>Host</b>	Rabbit
<b>Application</b>	IHC,
<b>Reactivity</b>	Human,Mouse,Rat

## Performance

<b>Conjugation</b>	Unconjugated
<b>Modification</b>	Unmodified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Form</b>	Liquid
<b>Storage</b>	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
<b>Buffer</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% New type preservative N.
<b>Purification</b>	Affinity purification

## Immunogen

<b>Gene Name</b>	PIK3R2
<b>Alternative Names</b>	p85; p85 beta; p85-BETA; P85B; P85B_HUMAN; Phosphatidylinositol 3 kinase; Phosphatidylinositol 3 kinase regulatory beta subunit; Phosphatidylinositol 3 kinase regulatory subunit beta; Phosphatidylinositol 3 kinase regulatory subunit polypeptide 2; Phosphatidylinositol 3 kinase, regulatory subunit, polypeptide 2 (p85 beta); Phosphatidylinositol 3-kinase 85 kDa regulatory subunit beta; phosphatidylinositol 3-kinase; Phosphatidylinositol 3-kinase regulatory beta subunit; Phosphatidylinositol 3-kinase regulatory subunit beta; Phosphoinositide 3 kinase regulatory subunit 2 (beta); Phosphoinositide 3 kinase regulatory subunit 2; Phosphoinositide 3 kinase regulatory subunit polypeptide 2 (p85 beta); Phosphoinositide 3 kinase regulatory subunit polypeptide 2; Phosphoinositide 3 kinase, regulatory subunit 2 (beta); Phosphoinositide 3 kinase, regulatory subunit 2 (p85 beta); PI3 kinase p85 beta subunit;

**Product Name: PI 3-Kinase p85 $\beta$  Rabbit Polyclonal Antibody**  
**Catalog #: APRab16105**

---

PI3 kinase p85 subunit beta; PI3-kinase regulatory subunit beta; PI3-kinase subunit p85-beta; PI3K; PI3K regulatory subunit beta; PIK3R 2; PIK3R2; polypeptide 2 (p85 beta); PtdIns 3 kinase p85 beta; PtdIns-3-kinase p85-beta; PtdIns-3-kinase regulatory subunit beta; PtdIns-3-kinase regulatory subunit p85-beta.

**Gene ID**

5296.0

**SwissProt ID**

O00459.Synthetic peptide from human protein at AA range: 1-60

## Application

**Dilution Ratio**

IHC-p 1:50-200, ELISA 1:10000-20000.

**Molecular Weight**

## Background


Phosphatidylinositol 3-kinase (PI3K) is a lipid kinase that phosphorylates phosphatidylinositol and similar compounds, creating second messengers important in growth signaling pathways. PI3K functions as a heterodimer of a regulatory and a catalytic subunit. The protein encoded by this gene is a regulatory component of PI3K. Two transcript variants, one protein coding and the other non-protein coding, have been found for this gene. [provided by RefSeq, Dec 2012],function: Binds to activated (phosphorylated) protein-tyrosine kinases, through its SH2 domain, and acts as an adapter, mediating the association of the p110 catalytic unit to the plasma membrane.,similarity: Belongs to the PI3K p85 subunit family.,similarity: Contains 1 Rho-GAP domain.,similarity: Contains 1 SH3 domain.,similarity: Contains 2 SH2 domains.,subunit: Heterodimer of a p110 (catalytic) and a p85 (regulatory) subunits.,

## Research Area


ErbB\_HER; Chemokine; Phosphatidylinositol signaling system; mTOR; Apoptosis\_Inhibition; Apoptosis\_Mitochondrial; Apoptosis\_Overview; VEGF; Focal adhesion; Toll\_Like; Jak\_STAT; Natural killer cell mediated cytotoxicity; T\_Cell\_Receptor; B\_Cell\_Antigen; Fc epsilon RI; Fc gamma R-mediated phagocytosis; Leukocyte transendothelial migration; Neurotrophin; Regulates Actin and Cytoskeleton; Insulin\_Receptor; Progesterone-mediated oocyte maturation; Type II diabetes mellitus; Aldosterone-regulated sodium reabsorption; Pathways in cancer; Colorectal cancer; Renal cell carcinoma; Pancreatic cancer; Endometrial cancer; Glioma; Prostate cancer; Melanoma; Chronic myeloid leukemia; Acute myeloid leukemia; Small cell lung cancer; Non-small cell lung cancer;

## Image Data


**Product Name: PI 3-Kinase p85 $\beta$  Rabbit Polyclonal Antibody**  
**Catalog #: APRab16105**



Immunohistochemical analysis of paraffin-embedded Human-testis, antibody was diluted at 1:100




Immunohistochemical analysis of paraffin-embedded Human-testis, antibody was diluted at 1:100



Immunohistochemical analysis of paraffin-embedded Human-colon, antibody was diluted at 1:100

**Product Name: PI 3-Kinase p85 $\beta$  Rabbit Polyclonal Antibody**  
**Catalog #: APRab16105**

---



Immunohistochemical analysis of paraffin-embedded Human-colon, antibody was diluted at 1:100

**Note**

For research use only.