

Product Name: Rab 20 Rabbit Polyclonal Antibody
Catalog #: APRab16744



Summary

Production Name	Rab 20 Rabbit Polyclonal Antibody
Description	Rabbit Polyclonal Antibody
Host	Rabbit
Application	WB,
Reactivity	Human,Rat,Mouse

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% New type preservative N.
Purification	Affinity purification

Immunogen

Gene Name	RAB20
Alternative Names	RAB20; Ras-related protein Rab-20
Gene ID	55647.0
SwissProt ID	Q9NX57.The antiserum was produced against synthesized peptide derived from human RAB20. AA range:115-164

Application

Dilution Ratio	WB 1:500 - 1:2000. ELISA: 1:10000..
Molecular Weight	28kD

Background


Product Name: Rab 20 Rabbit Polyclonal Antibody
Catalog #: APRab16744




function:Plays a role in apical endocytosis/recycling.,similarity:Belongs to the small GTPase superfamily. Rab family.,subcellular location:Highly enriched on apical endocytic structures in polarized epithelial cells of kidney proximal tubules.,tissue specificity:Low or absent expression in normal pancreas and stronger expression in 15 of 18 exocrine pancreatic adenocarcinomas (at protein level).,function:Plays a role in apical endocytosis/recycling.,similarity:Belongs to the small GTPase superfamily. Rab family.,subcellular location:Highly enriched on apical endocytic structures in polarized epithelial cells of kidney proximal tubules.,tissue specificity:Low or absent expression in normal pancreas and stronger expression in 15 of 18 exocrine pancreatic adenocarcinomas (at protein level).,

Research Area

Image Data




Western blot analysis of lysates from HeLa, MCF-7, HUVEC, and A549 cells, using RAB20 Antibody. The lane on the right is blocked with the synthesized peptide.




Western blot analysis of the lysates from HeLa cells using RAB20 antibody.

Product Name: Rab 20 Rabbit Polyclonal Antibody
Catalog #: APRab16744



Western Blot analysis of various cells using Rab 20 Polyclonal Antibody diluted at 1: 2000



Western Blot analysis of A549 cells using Rab 20 Polyclonal Antibody diluted at 1: 2000

Note

For research use only.