

Product Name: Ribosomal Protein L10L Rabbit Polyclonal Antibody
Catalog #: APRab17144

Summary

Production Name	Ribosomal Protein L10L Rabbit Polyclonal Antibody
Description	Rabbit Polyclonal Antibody
Host	Rabbit
Application	WB,ELISA
Reactivity	Human,Monkey

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% New type preservative N.
Purification	Affinity purification

Immunogen

Gene Name	RPL10L
Alternative Names	RPL10L; 60S ribosomal protein L10-like
Gene ID	140801.0
SwissProt ID	Q96L21.The antiserum was produced against synthesized peptide derived from human RPL10L. AA range:151-200

Application

Dilution Ratio	WB 1:500 - 1:2000. ELISA: 1:10000.
Molecular Weight	20kD

Product Name: Ribosomal Protein L10L Rabbit Polyclonal Antibody
Catalog #: APRab17144


Background

This gene encodes a protein sharing sequence similarity with ribosomal protein L10. It is not currently known whether the encoded protein is a functional ribosomal protein or whether it has evolved a function that is independent of the ribosome. This gene is intronless. [provided by RefSeq, Jul 2008],function:May play a role in compensating for the inactivated X-linked gene during spermatogenesis.,miscellaneous:This gene has no introns in its coding regions, and therefore was most likely produced by retrotransposition of the original X-linked gene during evolution.,similarity:Belongs to the ribosomal protein L10e family.,tissue specificity:Testis specific.,


Research Area

Ribosome;

Image Data




Western blot analysis of lysates from Jurkat, COS7, and HeLa cells, using RPL10L Antibody. The lane on the right is blocked with the synthesized peptide.



Western blot analysis of the lysates from HUVEC cells using RPL10L antibody.

**Product Name: Ribosomal Protein L10L Rabbit
Polyclonal Antibody
Catalog #: APRab17144**



Western Blot analysis of various cells using Ribosomal Protein L10L Polyclonal Antibody diluted at 1 : 1000

Note

For research use only.