

Product Name: Ribosomal Protein S13 Rabbit Polyclonal Antibody
Catalog #: APRab17174

Summary

Production Name	Ribosomal Protein S13 Rabbit Polyclonal Antibody
Description	Rabbit Polyclonal Antibody
Host	Rabbit
Application	IHC,ELISA
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% New type preservative N.
Purification	Affinity purification

Immunogen

Gene Name	RPS13
Alternative Names	RPS13; 40S ribosomal protein S13
Gene ID	6207.0
SwissProt ID	P62277.The antiserum was produced against synthesized peptide derived from human RPS13. AA range:71-120

Application

Dilution Ratio	IHC 1:100-1:300 ELISA: 1:40000
Molecular Weight	17kD

Background


Product Name: Ribosomal Protein S13 Rabbit Polyclonal Antibody
Catalog #: APRab17174

Ribosomes, the organelles that catalyze protein synthesis, consist of a small 40S subunit and a large 60S subunit. Together these subunits are composed of 4 RNA species and approximately 80 structurally distinct proteins. This gene encodes a ribosomal protein that is a component of the 40S subunit. The protein belongs to the S15P family of ribosomal proteins. It is located in the cytoplasm. The protein has been shown to bind to the 5.8S rRNA in rat. The gene product of the E. coli ortholog (ribosomal protein S15) functions at early steps in ribosome assembly. This gene is co-transcribed with two U14 small nucleolar RNA genes, which are located in its third and fifth introns. As is typical for genes encoding ribosomal proteins, there are multiple processed pseudogenes of this gene dispersed through the genome. [provided by RefSeq, Jul 2008],similarity:Belongs to the ribosomal protein S15P family.,

Research Area

Ribosome;

Image Data



Immunohistochemistry analysis of paraffin-embedded human breast carcinoma tissue, using RPS13 Antibody. The picture on the right is blocked with the synthesized peptide.

Note

For research use only.